# **Contec® DualClean™ Quiltec® Cover**

# Sealed edge polyester cover for use with DualClean Softwall Strip Cleaner Frame

The DualClean Softwall Strip Cleaner is a unique system that enables one operator to easily and quickly clean both sides of 4"-12" (10 - 30cm) plastic strip curtains in a single pass there by reducing time to clean softwall strips. The system consists of a cover that is installed over a stainless steel frame with two curved-surface arms.

Made of 2-ply knitted polyester, the Quiltec Cover for the DualClean Softwall Strip Cleaner is compatible with a wide variety of solvents and disinfectants. It is ideally suited for more critical cleanroom environments.

The DualClean Softwall Ouiltec Cover is laundered in an ISO Class 4 Cleanroom.

For best performance, use saturated with suitable solvent or disinfectant.

Covers are fully autoclavable with a recommended cycle of 121°C for 30 minutes.

# Features & Benefits

• Extremely clean and durable

**Performance Characteristics** 

- Can be sterilized by most common methods making it ideal for critical environments
- Improves cleaning performance by trapping and removing contaminants
- Compatible with most common solvents and disinfectants

















### **Product Information** Cover material 100% knitted polyester Construction 2-layer, ultrasonically bonded, interlock knit

	, , , , , , , , , , , , , , , , , , , ,
ISO Class cGMP Grade	4-8 A/B - for sterile version, C/D - for non-sterile version
Compatible Mop Frames	2778
Shelf Life	For Sterile Products: 2 years from manufacturing date

See other side for additional information.

Copyright © 2019 Contec, Inc. All rights reserved.

062519

# **Technical Data**

Attribute (units)	Typical Value*	Test Method**	
Basis weight; nominal (g/cover)	31	Contec Method	
Sorbent capacity; (mL/cover)	92.9	CTM004-03	

#### \*Notes

- a) The data shown are typical values and should not be used as product specifications.
- b) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- c) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

#### \*\*Test Methods:

- 1) CTM Refers to Contec Test Method
- 2) IEST-RP-CC004.3, Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments, Institute of environmental Sciences and Technology, Rolling Meadows IL.

# **Products**

Part Number	Description	EA/OB	OB/CS	EA/CS
DCQT0001	Contec DualClean Quiltec Cover, knit polyester, 26.3" x 2.77" (67 x 7cm)	8	10	80

Part Number	Description	EA/OB1	OB1/OB2	OB2/OB3	OB3/CS	EA/CS
DCQT0002	Contec DualClean Quiltec Cover, knit polyester, 26.3" x 2.77" (67 x 7cm) sterile, triple-bagged	4	1	1	20	80

## **Additional Information**

For Sterile products, Certificates of Conformance and Sterility are available at www.contecinc.com/Resources/Certificate Search Products are Validated Sterile by Gamma irradiation to a Sterility Assurance Level of 10<sup>-6</sup> according to AAMI Guidelines Contec holds ISO 9001:2015 registration. Additional Product information available at www.contecinc.com.

